



23 CLINICAL EVALUATORS

215 TOTAL USES

88% CLINICAL RATING

Key features: Fast-Setting VPS Impression Material • High Detail
• Color Contrast

Description

PERFIT is a polyvinylsiloxane impression material which is:

- Available in 4 viscosities: Heavy, Regular, Light, and Putty
- Fast setting: sets in 3 minutes intraorally

Indication

PERFIT impression materials are used to make final impressions for all indirect restorative procedures, including:

- crown
- bridge
- inlays and onlays
- implantology
- orthodontics
- removable appliances



Clinical Tip

- Keep field as dry as possible. Be sure light body covers tooth with no air bubbles or gaps.
- The light body works really well for a denture reline impression material.

"GREAT DETAIL AND HANDLING. EXCELLENT HYDROPHILIC PROPERTY."

Evaluators' Comments

"The ease of loading and consistency of this impression material are very good."

"No voids, light body is more runny and goes easy to the required area to be captured."

"The color of the light body made it easy to read margin details."


"I loved the set-up time!"

"I prefer a brighter color of light body."


"Even though this was fast set, I would like a faster set."

"Label cartridges with set time."

Evaluation Summary:



Compared to Competitive Products:



Consultants who would:
78% Recommend to a colleague

Consultants who would
want to stock in office:
43% Yes, instead of current product