



35 CLINICAL EVALUATORS

786 TOTAL USES

91% CLINICAL RATING

Key features: Autoclavable • Multi-Use • Color-Coded



Description

DIA-BURS are:

- Autoclavable, reusable diamond burs available in five grits
- Available in numerous shapes and sizes

Indications

- All operative procedures for direct and indirect restorations requiring reduction of tooth structure or restorative materials

Unique Attributes

- Diamonds plated with electroplating technology; this minimizes or eliminates diamond particles separating from the shank after repeated use
- Various shapes and sizes of diamond particles enhance sharpness



Clinical Tip

- Although I don't use the shape regularly, I found the SF-41 diamond very helpful in several applications including crown prep and slots for extra coronal splinting

"MAINTAINS SHARPNESS, EVEN WHEN REMOVING FAILED ZIRCONIA RESTORATIONS."

Evaluators' Comments

"Consistency of the cut stayed the same, even after 3-4 uses."

"Great, efficient cutting diamonds. They cut well even after autoclaving."

"Diamond particles did not fall off like other brands."

"I was very impressed at how sharp they stayed."

"Packaging did not hold burs well. Recommend a change to single packaging."

Consultants who would:

91%

Recommend to a colleague

Consultants who would want to stock in office:

15%

Yes, instead of current product


65%

Yes, in addition to current product

12%

I might want to order this product for certain cases

Evaluation Summary:



Compared to Competitive Products:

