



25 CLINICAL EVALUATORS

659 TOTAL USES

91% CLINICAL RATING

Key features: Universal adhesive • Two-bottle system

Description

G-2 Bond Universal is a universal bonding agent with:

- Choice of preferred etching mode: self-adhesive, selective etch, or total etch
- High bond strength and durability
- Multiple indications
- Resistance to marginal staining
- Easy dispensing with bottle covers that control the pressure for precise drops
- Low technique sensitivity
- HEMA free-formulation

Indications

- Direct bonding of light-cured composites and compomers to tooth structure
- Bonding of dual-cured core build-up composites to tooth structure as long as these materials are light-cured
- Intraoral repair of composite, metal-based and zirconia/alumina-based restorations
- Intraoral repair of ceramic restorations in combination with a silane coupling agent
- Treatment of hypersensitive teeth
- Sealing of tooth preparation (cavity or abutment) for indirect restorations
- Treatment of exposed root surfaces
- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges, and veneers
- Post-cementation use



Unique Attributes

- This two-bottle system helps to reduce technique sensitivity that can lead to post-op sensitivity and restoration failures. Evaluator comments confirmed that using this material reduced the incidence of post-op sensitivity.
- The precise dropper system reduces wasted material.

Clinical Tip

- Protect from ambient light when doing multiple restorations, as it sets up quickly in the well.

"GREAT FOR REDUCING SENSITIVITY."

Evaluators' Comments

"A universal like **G2 Bond Universal** should be in every office."

"Great for reducing sensitivity."

"I love the precise dispensing."

"Worked great to repair porcelain on crowns."

"Good film thickness and flowability."

"I typically use a desensitizer to minimize post op sensitivity. I did not in this case and did not have any patients report sensitivity."

"Cured quickly in the well and I had to use multiple drops when placing several restorations."

"Great product but would have preferred if it were 1 step."

Consultants who would:

76% Recommend to a colleague


Consultants who would want to stock in office:

16% Yes, instead of current product

44% Yes, in addition to current product

12% I might want to order this product for certain cases

Evaluation Summary:



Compared to Competitive Products:

