



31 CLINICAL EVALUATORS

218 TOTAL USES

93% CLINICAL RATING

Key features: Esthetic fiber post • Available in six sizes • Pre-primed
• Matching drills

Description

Super Fiber Posts®:

- "Total internal reflection" provides excellent light transmission for effective curing
- Elasticity similar to dentin with low modulus of elasticity
- High radiopacity
- Color-coded stopper indicating size of post
- Pre-primed surface
- Packaged individually in blister packs

Indications

- One-step post and core buildup procedures
- Use with light-cured or dual-cured core-buildup material

Clinical Tips

- Use other post drills to just smaller than post you will use, then use coordinating drill from kit.
- Remember to remove 2 mm of the coronal dentin per the instructions.
- On some occasions I removed the stoppers as they were in the way.

Unique Attributes

- Pre-primed surface
- Hygienic packaged in blister packs



"THE
RADIO-PACITY
WAS THE BEST
I'VE SEEN!"

Evaluators' Comments

"Love the ease of use and individual packs, it makes placement quick and easy."

"Cutting excess was super easy."

"Cementation was easy since the posts are pre-primed."

"The post prep drills were easy to use and left a clean surface."

"The posts adapted very well to the canals and were very radiopaque."

"I liked the posts a lot but was disappointed in the drills."

"The coronal portion of the post was substantial in size compared to others."

"I liked the blister packs and the color coding."


"Great light transmission and the instructions for use are good."

"The drills did not cut dentin very well."


"The packaging is not very convenient or sturdy."

"I would like longer posts."

Evaluation Summary:



Compared to Competitive Products:



Percentage of Consultants Who Would:

58% Recommend in addition to product they currently use

32% Recommend instead of product they currently use

10% Not recommend